



Recommendations of SEPAR

Pediatric Bronchoscopy Guidelines[☆]

Normativa de broncoscopia pediátrica

Javier Pérez-Frías^{a,*}, Antonio Moreno Galdó^b, Estela Pérez Ruiz^a,
María Isabel Barrio Gómez De Agüero^c, Amparo Escribano Montaner^d, Pilar Caro Aguilera^a

^a Unidad de Neumología Infantil, Hospital Universitario Materno Infantil Carlos Haya, Facultad de Medicina, Málaga, Spain

^b Unidad de Neumología Pediátrica y Fibrosis Quística, Hospital Universitari Vall d'Hebron, Universitat Autònoma de Barcelona, Barcelona, Spain

^c Sección de Neumología Pediátrica, Hospital Infantil La Paz, Madrid, Spain

^d Unidad Neumología Infantil, Hospital Clínico Universitario, Valencia, Spain

ARTICLE INFO

Article history:

Received 18 March 2011

Accepted 2 April 2011

Available online 8 September 2011

Keywords:

Lexible bronchoscopy

Children

Bronchial biopsy

Transbronchial biopsy

Bronchoalveolar lavage

Rigid bronchoscopy

Palabras clave:

Broncoscopia flexible

Niños

Biopsia bronquial

Biopsia transbronquial

Lavado broncoalveolar

Broncoscopia rígida

Introduction

The development that bronchoscopy has experienced in the recent years has been spectacular, especially in the field of pediatrics. Both the diagnostic and the therapeutic applications of bronchoscopy have increased considerably. This expansion has been mainly based on the technological developments in various areas, such as instrumentation, optics, fiber optics, light sources, video-electronics and anesthesia techniques, to name just an important few.

The aim of these guidelines is to improve, facilitate and unify the process of bronchoscopy in pediatrics. The recommendations have been established in keeping with the GRADE system¹ (Table 1). In

the aspects in which the scientific evidence is insufficient, we have included the recommendations agreed upon by consensus of the authors.

Flexible Bronchoscopy

Indications

Flexible bronchoscopy (FB) enables us to obtain anatomical and dynamic information of the airways and to perform cytological and microbiological studies.² Its indications arise with the need to respond to symptoms or radiological anomalies that cannot be explained by non-invasive methods³ or to obtain samples from the lower airways⁴ (Table 2). In children, compared with rigid bronchoscopy (RB), it is not a good therapeutic tool,⁵ although it does facilitate procedures that can resolve some of the problems in children's airways.⁶

[☆] Please cite this article as: Pérez-Frías J, et al. Normativa de broncoscopia pediátrica. Arch Bronconeumol. 2011;47:350–60.

* Corresponding author.

E-mail address: fjperez@uma.es (J. Pérez-Frías).

Table 1
Classification of the Recommendations and Quality of the Evidence According to the GRADE System.^a

Grade of Recommendation	Quality of Supporting Evidence	Implications
Strong recommendation ^b High-quality evidence	Well-performed RCT or exceptionally well-done, unbiased OS	Can apply to most patients under most circumstances
Strong recommendation ^b Moderate-quality evidence	RCT with important limitations or unusually strong unbiased OS	Can apply to most patients under most circumstances
Strong recommendation ^b Low-quality evidence	Evidence for at least one critical outcome from OS or from RCT with serious flaws or indirect evidence	May change when higher quality evidence becomes available
Strong recommendation ^b Very-low-quality evidence	Evidence for at least one of the critical outcomes from unsystematic clinical observations or very indirect evidence	May change when higher quality evidence becomes available
Weak recommendation ^c High-quality evidence	Well-performed RCT or exceptionally strong, unbiased OS	May differ depending on circumstances or patients
Weak recommendation ^c Moderate-quality evidence	RCT with important limitations or unusually strong unbiased OS	Alternative approaches likely to be better for some patients under some circumstances
Weak recommendation ^d Low-quality evidence	Evidence for at least one critical outcome from OS or from RCT with serious flaws or indirect evidence	Other alternatives may be equally reasonable
Weak recommendation ^e Very-low-quality evidence	Evidence for at least one of the critical outcomes from unsystematic clinical observations or very indirect evidence	Other alternatives may be equally reasonable

RCT: randomized control trials; OS: observational studies.

^a Adapted by Schünemann et al.¹

^b The benefits clearly outweigh the drawbacks or vice versa.

^c The benefits are balanced with the drawbacks.

^d Uncertainty in the estimation of the benefits or drawbacks. The benefits may be balanced with the drawbacks.

^e Major uncertainty in the estimation of the benefits or drawbacks, the benefits may or may not be balanced with the drawbacks.

Table 2
Indications for Flexible Bronchoscopy in Pediatrics.

1. Exploration of the airways
<i>Persistent stridor</i>
<i>Persistent wheezing</i>
<i>Hemoptysis</i>
<i>Suspicion of a foreign body</i>
<i>Persistent or recurring atelectasis</i>
<i>Persistent or recurring pneumonia</i>
<i>Localized pulmonary hyperlucency</i>
<i>Problems related with the artificial airways</i>
<i>Miscellaneous (large burns, phonatory anomalies)</i>
2. Obtaining biological samples (bronchoalveolar lavage, bronchial biopsy, bronchial brushing)
<i>Pneumonia in immunosuppressed patients</i>
<i>Chronic interstitial pneumonia</i>
<i>Pneumonitis due to hypersensitivity</i>
<i>Pulmonary hemosiderosis</i>
<i>Eosinophilic pneumonia</i>
<i>Other (sarcoidosis, alveolar proteinosis, histiocytosis)</i>
<i>Endoluminal obstructive pathology</i>
<i>Aspiration lung syndromes</i>
3. Therapeutic applications
<i>Aspiration of endobronchial secretions</i>
<i>Instillation of medication</i>
<i>Difficult or selective intubations</i>
<i>Management of the foreign body combined with rigid bronchoscope</i>

Examination of the Airways

Persistent Stridor

Its main cause is laryngomalacia,⁷ which tends to disappear within the first year and therefore does not usually require endoscopic revision, except in the cases of parental anxiety. In the cases of atypical presentation, biphasic character, prolonged persistence, suffocation crisis, eating difficulties, association with syndromes or malformations, history of intubation or diagnosis of severe laryngitis in children under the age of 6 months, a complete exploration is recommended⁸ as there may be anatomical, congenital or acquired anomalies⁹ (strong recommendation; high-quality evidence).

Suspicion of a Foreign Body

The existence of respiratory symptoms and/or recurrent or persistent radiological images can be associated with the presence of an unnoticed foreign body. Under these circumstances, RB has a high rate of false negatives. In contrast, FB has a diagnostic exactitude of 100%. Therefore, if a diagnosis is mistaken, it is recommended to first perform FB¹⁰ (strong recommendation; high-quality evidence). Once extracted, FB can be very useful for confirming the absence of residues or of fragmented material, as occurs in the case of nuts.

Difficult-to-Control, Persistent Wheezing

FB is recommended in cases of difficult-to-control wheezing bronchitis, especially in small children, associated with asymmetries upon auscultation and/or radiological alterations¹¹ as they could be due to a foreign body or structural anomalies of the tracheobronchial wall¹² (strong recommendation; high-quality evidence).

Hemoptysis

Hemoptysis is an uncommon and alarming symptom. In the absence of a known explanation (otorhinolaryngological infections or bronchiectasis), FB is required¹¹ (strong recommendation; high-quality evidence). The most frequent causes are problems related with artificial airways, tracheostomies, primary alveolar hemorrhage, congenital vascular anomalies and infectious or inflammatory endobronchial pathologies, while tumors are the exception.¹³

Persistent/Recurring Atelectasis

Persistence of more than 6 weeks, together with unexplained symptoms, makes FB recommendable¹⁴ (weak recommendation; high-quality evidence). The most frequent findings are mucus plugging, foreign bodies, extrinsic compression in cases of congenital cardiopathy,^{15,16} inflammatory granulation tissue or endobronchial tuberculosis.^{17,18}

Persistent/Recurring Pneumonia

Middle lobe syndrome is a relatively frequent entity in pediatrics due to the anatomical characteristics of this bronchus. The

observation of bacterial growth, a predominance of neutrophils in the bronchoalveolar lavage and evolution towards bronchiectasis in more than 50% of children with non-atopic wheezing¹⁹ has renewed interest in this syndrome.²⁰ It has been proposed that the syndrome, if it cannot be resolved with treatment, indicates a study with FB²¹ (strong recommendation; high-quality evidence).

Localized Hyperlucency

When its presence is not associated with congenital or post-infectious causes, localized hyperlucency requires carrying out FB as there may be air trapping due to a valvular mechanism, which is usually secondary to an intrinsic obstruction (foreign body, bronchomalacia,²² granulation tissue) or extrinsic compression due to adenopathies, vessels that are aberrant, anomalous and/or increased in size, or congenital or acquired mediastinal masses²³ (strong recommendation; high-quality evidence).

Problems Related With Artificial Airways²⁴

FB is essential for the diagnosis of problems that arise during in/extubation and for monitoring and following patients with tracheostomy²⁵ (strong recommendation; high-quality evidence). The most frequent findings are laryngeal edema, subglottic stenosis and tracheobronchial granulation tissue secondary to trauma caused by endotracheal tubes and cannulae or by repeated aspirations.²⁶

Miscellaneous²⁷

Phonatory abnormalities: abnormal or persistent, aphonic crying can be due to the unilateral paralysis of a vocal cord generally associated with trauma during birth, hydrocephalus or post-operative cardiopathy, juvenile papillomatosis or aspiration during birth. The association with stridor warrants carrying out FB with a complete examination of the airway²⁸ (strong recommendation; high-quality evidence). Other situations indicating FB include: persistent cough without diagnosis by non-invasive techniques (as it may be due to entities such as bacterial malacia or bronchitis.)²⁹ and the evaluation of injury secondary to burns or inhalation of toxic substances.³⁰

Bronchoalveolar Lavage and Biopsy Samples

FB is widely used in researching lung infiltrates that are unclarified, localized or diffuse, with interstitial, alveolar, miliary or nodular pattern, when diagnosis is not possible by means of other less-invasive methods (strong recommendation; high-quality evidence).

Therapeutic Procedures

Aspiration of Secretions

FB can be useful for resolving atelectasis due to the retention of secretions or mucus plugging. The percentage of radiological re-expansion is variable, generally no higher than 50%. In cystic fibrosis³¹ (weak recommendation; low-quality evidence) or in plastic bronchitis (moderate recommendation; low-quality evidence) endoscopic lavage and aspiration can be indicated together with the instillation of mucolytics³² (MESNA or DNasa).³³ In alveolar proteinosis, repeated BAL is the recommended treatment³⁴ (strong recommendation; high-quality evidence).

Difficult and Selective Indications

Bronchoscopy can act as a guideline for intubation in cases of craniofacial anomalies and polymalformative syndromes,³⁵ or in selective bronchial intubation (strong recommendation; high-quality evidence).

Extraction of Foreign Bodies

Although there are publications that endorse the good performance of FB,^{36,37} it is usually a complicated procedure in small children. In these patients, RB is preferred as it offers the advantages of general anesthesia, assisted ventilation, larger instruments and a greater variety of accessories³⁸ (strong recommendation; high-quality evidence). The ideal procedure would be to initiate with FB, which allows for greater reach in the exploration and the identification of the foreign body, extraction with RB, and a final revision with FB in order to rule out a residual foreign body.¹⁰ In some cases in which the clinical and/or radiological information is conclusive, the procedure may be directly initiated with RB.¹⁰

Other Procedures

FB can be of use in the peri-operative management of tracheoesophageal fistulas,³⁹ both isolated or recurrent, as it allows for its identification and canalization, facilitating the surgical approach⁴⁰ (weak recommendation, low-quality evidence). In other techniques such as bronchial permeabilization or the placement of prosthesis, the small size of the FB can be either a complication or an impediment, and therefore RB should be used⁴¹ (strong recommendation; high-quality evidence).

Contraindications

The indication for FB should be individualized, evaluating the risk/benefit for each patient. Some absolute contraindications impede performing bronchoscopy: severe refractory hypoxemia, hemodynamic instability, uncorrected hemorrhagic diathesis or the lack of authorization for the procedure by the parent or guardian. There are some relative contraindications determined by the experience of the team or the level of critical care of the hospital: very premature newborns, congenital cyanotic cardiopathies with an increase in bronchial collateral circulation, severe pulmonary hypertension or coagulation alterations.

Complications

Pediatric FB is in general a safe procedure. The possible complications depend on the patient's risk factors, the type of procedure done, inadequate sedation/anesthesia, choice of instruments of an improper size and the inexperience of the bronchoscopist. The main complications are as follows.

Nasal Trauma and Epistaxis

These are the most frequent. Compression and/or topical adrenalin are sufficient for its control.

Desaturation and Hypoxemia

These may be due to causes inherent of the procedure and patient, or to the excess or deficit of sedation. FB increases the resistance of the airway which, together with sedation, can motivate hypoxemia and hypercapnia, which can persist 15–20 min after the exploration. BAL and prolonged aspiration maneuvers, mainly in newborns and infants, favor alveolar collapse and increase its risk. Due to vagal irritation, or in children with bronchial hyperreactivity, there may also be some degree of bronchospasm. These changes can be attenuated by maintaining a stable hemodynamic situation, choosing instruments of the proper size, previously administering bronchodilators and supplemental oxygen, optimizing the sedation-anesthesia and performing the procedure quickly.⁴²

Cough and/or Bronchospasm

These may be related with not enough sedation or the instillation of topical anesthesia.

Table 3
Pediatric Flexible Bronchoscopes and Video-Assisted Endoscopes.

External diameter, mm	WC, mm	BF	PC/PCB	ETT, mm	Utility
4.9 (FB)	2.2	+	+	6	>7–8 years or >20 kg
4.9 (VE)	2.0				
3.6 (BB)	1.2	+	–	5	Standard pediatric
3.8 (VE)					
2.8 (FB and VE)	1.2	+	–	4	Newborns–infants–children
2.2 (FB)	–	–	–	3	Newborns and infants <6 months

WC: diameter of the working channel; PC/PCB: protected catheters and protected catheter brushes; FB: flexible bronchoscope; BF: biopsy forceps; ETT: minimal internal diameters of endotracheal tubes for the passage of the bronchoscope; VE: video endoscopes.

Trauma and Obstruction of the Airway

Also related with the choice of bronchoscope of improper size or forcing the way through a stenotic area, causing edema, hypersecretion and greater obstruction.

Hemorrhage

Aspiration and brushing or biopsies usually cause small transitory hemorrhages facilitated by the existence of granulation tissue, tumors, bronchiectasis or hemorrhagic diathesis. Severe hemorrhage is exceptional, and can be associated with overly vascularized endobronchial lesions, which are rare in children, BAL in patients with coagulation disorders and especially transbronchial biopsy.

In general, hemorrhaging will let up spontaneously in a few minutes, but in cases of severe hemorrhage, lavage with cold physiological saline solution can be of use (5 ml aliquots) or with adrenalin at 1/20 000–1 ml of adrenalin at 1/1000, diluted in 20 ml of physiological saline solution (strong recommendation; low-quality evidence). If the hemorrhage persists, it is recommended to place the patient in lateral decubitus, with the bleeding lung in a declining position, and to insert the bronchoscope in a segmental or subsegmental bronchus, aspirating continuously for 3–5 min in order to produce a collapse of the distal bronchial walls (strong recommendation; low-quality evidence). On rare occasions, compression is needed with a 3 Fr Fogarty catheter or selective guided intubation. In these cases, RB could be necessary to explore the area of the bleeding and to aspirate the large clots. In one study, it has been suggested that the endobronchial instillation of procoagulation substances, such as tranexamic acid, can be effective, although the evidence for its recommendation is still limited⁴³ (weak recommendation; low-quality evidence).

Pneumothorax

Pneumothorax is a rare complication although it can occur with transbronchial biopsy. In other cases, the misuse of the bronchoscope (uncontrolled pressure with the tip, abrupt movements, especially if there is coughing, or using either O₂ or pressured air through the working channel) could cause them.⁴⁴

Fever and Infections

Fever is frequent, appearing in up to 15% of the procedures, especially if BAL had been carried out. It is related with the release of cytokines or with transitory bacteriemia, mainly in immunosuppressed patients. More rarely, FB could give rise to infectious complications, both for the patient (dissemination to healthy areas of the lung or transmission of contaminated material from a previous exploration) as well as for the health-care professional (transmission of *M. tuberculosis*).

Instruments

In pediatrics, different models of bronchoscopes are used depending on the age of the child or the size of the endotracheal tube through which it is introduced (Table 3). The external diame-

Table 4
Biopsy Forceps Available for Pediatric Bronchoscopes.

Diameter, mm	Characteristics
1 ^a	Round, smooth edge Oval, smooth, disposable Oval, alligator
1.2 ^b	Oval spoon, mouse-tooth
1.5 ^a	Oval, smooth, disposable Dentate spoon
1.8 ^a	Round, smooth edge Oval, smooth, fenestrated, disposable Oval, smooth, fenestrated, disposable, with a needle Dentate spoon
2 ^b	<i>Re-usable</i> Fenestrated, mouse-toothed alligator Round, fenestrated with/without needle Oval, fenestrated Oval, fenestrated alligator <i>Disposable</i> Oval, fenestrated Oval, fenestrated with needle Oval, fenestrated alligator Oval, fenestrated alligator with needle

^a Refers to the diameter of the forceps.

^b Refers to the diameter of the working channel.

ter of the flexible tube and the working channel vary according to the type of bronchoscope used. There are four different bronchoscopes used in children, sized according to their external diameter: 2.2 mm, 2.8 mm, 3.5 mm and 4.9 mm.⁴⁵ The smallest (2.2 mm) does not have a working channel and is therefore only used for inspecting the airway.⁴⁶ The 2.8 mm and 3.5 mm scopes have a 1.2 mm channel. The 4.9 mm scope has a 2.2 mm channel, which enables protected brushing and biopsies and is more effective in aspirating secretions. This model can be used in children either over the age of 7 or weighing more than 20 kg. Currently, bronchoscopy is being substituted with video-assisted bronchoscopy, which provides greater resolution and storage of the images in digital format.

Accessories

- Biopsy forceps: the most often used are alligator forceps and they are used in bronchial as well as transbronchial biopsies. There are other models with smooth edges, with a central needle, fenestrated or not, revolving or swinging. All these are shown in Table 4.
- Brushes: used for bronchial cytology studies and to test signs of inflammation, infection and malignancy.
- Transbronchial aspiration needles: used for diagnosing lung nodules or lesions near the bronchial tree. They are not frequently used in children.
- Accessories for extracting foreign bodies: although the standard norm is not to extract foreign bodies with FB, it is useful to have the equipment available in order to recover either small fragments or those from more distal areas that are not accessible with

Table 5
Accessories for the Extraction of Foreign Bodies With Flexible Bronchoscope.

Diameter, mm	Characteristics
1 ^a	Helical basket, 4 wires
1.2 ^b	Helical baskets, 3 and 4 wires 3-pronged forceps Wire loop
1.3 ^a	Basket, 6 wires
1.5 ^a	Handle
1.8 ^a	Helical basket, 4 wires Mouse-toothed forceps V forceps 3-pronged forceps
2 ^b	Basket, 4 wires V forceps Mouse-toothed forceps Rubber forceps

^a Refers to the diameter of the forceps.

^b Refers to the diameter of the working channel.

the rigid bronchoscope. There are 3 models: mini-forceps, forceps and extraction baskets (Table 5).

- Balloon catheter: used to plug hemorrhagic lesions.
- Endotracheal tube adapter: used in intubated patients to allow for simultaneous ventilation.

Description of the Procedure

If possible, FB should be performed in patients under mild sedation, enough to tolerate the technique and the disorientation, permitting spontaneous breathing and thus reducing the risk of hypoxia. The previous conditions of the child will determine the type of access and the characteristics of the exploration.

Prior to carrying out the technique, the procedure will be explained, in comprehensible terms, to the child and/or parents and informed consent will be obtained. This will reduce the anxiety of family members and increase collaboration. A detailed clinical history should be available, including any circumstances that could affect the exploration (bronchial hyperresponsiveness, obstructive sleep apnea syndrome, allergies, cardiopathy and coagulation alterations).

The procedure is simple, well-tolerated and generally does not require hospitalization. The child should not ingest water during the 2 h previous to the technique nor any food during 6 h.

Options for Introducing the Flexible Bronchoscope

- Nasal: this requires IV sedation and topical anesthesia of the nostrils (2% lidocaine).
- Nasal with a facial mask: provides O₂ at 100% and continuous positive airway pressure support.
- Oral through a laryngeal mask: enables exploration of the lower airways from the glottic region.
- Endotracheal tube or tracheostomy: only allows for the observation of the lower airway.
- Rigid bronchoscope: its larger diameter accommodates the passage of the flexible bronchoscope.

Procedure

The operator sits at the head of the patient, with the gurney in low position in order for the disposition of the equipment to avoid a curve. Forced curves in the scope can damage its fibers and make it difficult to handle. If the access is oral, a mouth prop or bite block is placed in the patient's mouth, which is not necessary if the access is through the nasopharynx. The bronchoscope is chosen to be proper for the child's age, and the distal end is lubricated. It is connected to

the aspiration system, with pressure between 25 and 120 cmH₂O, interspersing appropriate receptors for sample collection.

If the access is nasal, the anatomy and the functionality of the pharyngeal and laryngeal structures (sublingual glands, tonsils, arytenoids, epiglottis and vocal cords) are studied. In cases where the larynx cannot be visualized (flaccidity, secretions), the line of air or oxygen can be connected to the working channel, at 1–2 l/min, in order to exert positive pressure to clean or open the area. It is not advisable to conduct such a procedure placing the bronchoscope below the larynx given the risk for pneumothorax (weak recommendation; low-quality evidence).

The translaryngeal passage is done by centering the end of the bronchoscope in the angle of the anterior corner of the vocal cords, introducing it by means of posterior flexion when the patient inhales. In order to make the passage easier and to prevent the appearance of laryngeal spasm, a local dose of 1% lidocaine can be administered through the working channel. After reaching the subglottic space, new doses of lidocaine can be instilled in the lower airways, without ever surpassing the calculated maximal dose. After each instillation, a small volume of air should be introduced through the bronchoscope channel to completely empty it.

During the procedure, it is necessary to observe the axis of the trachea and its movement during breathing, the presence of cartilage, *pars membranacea*, carina, possible compression-displacement, pulse areas, etc. The exploration of the bronchial tree should be systematic, complete and sequential.⁴⁷ In newborns weighing less than 1 kg, the scope may be introduced into the main bronchi and the opening of the lobes may be observed. In newborns weighing more than 2.5 kg, practically all the segments can be examined, except those of the upper lobes. At more than 6–8 kg, there should be no problems to explore all the segments. The airway should be checked for possible anatomical anomalies, while observing the appearance of the mucosa (pale, erythematous, friable, thin or thick), the characteristics of the secretions (scant, moderate or abundant; localized or diffuse; mucous, purulent, hemorrhagic), and proceed with the taking of samples (secretions or tissues) necessary.

Post-Bronchoscopy Care

The patient should be placed in a comfortable position. The patient should be watched for the appearance of complications. Full consciousness needs to be recuperated and oral tolerance should be checked before making the decision to release the patient. No food or drink should be offered until the swallowing reflex reappears (2 h after the procedure).

Bronchoalveolar Lavage

BAL allows for the evaluation of the cell and biochemical constituents of the epithelial surface of the lower respiratory tract or the presence of microorganisms, by means of the instillation and later aspiration of liquid in one or several lung segments.

Pre-requisites

Chest radiographs and, if possible, computed tomography in order to outline the most ideal segment for BAL. In certain patients, platelet count higher than 60 000/mm³ and Quick index >50%.

Technique

If bacteriological studies are going to be ordered, BAL should be done first, before exploring the bronchial tree, avoiding aspirating secretions with the aim of impeding, or reducing the contamination of the suction channel by the oropharyngeal flora (strong recommendation; low-quality evidence). Likewise, lavage should be done before other techniques (biopsy, brushing or needle aspiration) that could cause hemorrhages and affect the results. Refrain from the

Table 6
Differential Cell Count in the Bronchoalveolar Lavage of Different Studies in Normal Children.

	Clement et al. ⁴⁹	Ratjen et al. ⁵⁰	Riedler et al. ⁵¹	Midulla et al. ⁵²	Tessier et al. ⁵³
Number of children	11	48	18	16	11
Age, years	1–15	3–5	1–10	2–3	4–16
Number of aliquots	6	3	3	2	6
Volume of physiological saline	10% FRC	3 ml/kg	3 ml/kg	20 ml	10% FRC
Recovered BAL liquid, %					
Mean ± SD	Not done	58±15	Not done	43.1±12.2	69.7±9.6
Median	Not done	Not done	62.55	42.5	68
Range	Not done	Not done	42.5–71.5*	20–65	52–87
×10 ⁴ cells/ml					
Mean ± SD	25.5±13.6	10.3±11.1	Not done	59.9±32.9	35.1±18.4
Median	24	7.3	15.5	51	30.5
Range	7.0–50.0	0.5–57.1	7.5–25.8*	20–130	9–68
Macrophages, %					
Mean ± SD	89.7±5.2	81.2±12.7	Not done	86±7.8	89.9±5.5
Median	89	84	91	87	92.5
Range	85–97	34.6–94	84.2–94*	71–98	77–98
Lymphocytes, %					
Mean ± SD	8.7±4.6	16.2±12.4	Not done	8.7±5.8	8.9±5.6
Median	10	12.5	7.5	7	8
Range	1–17	2–61	4.7–12.8*	2–22	2–22
Neutrophils, %					
Mean ± SD	1.3±0.9	1.9±2.9	Not done	5.5±4.8	1.2±1.2
Median	1	0.9	1.7	3.5	1
Range	0–3	0–17	0.6–3.5*	0–17	0–3
Eosinophils, %					
Mean ± SD	Not done	0.4±0.6	Not done	0.2±0.3	0
Median	Not done	0.2	0.2	0	0
Range	Not done	0–3.6	0–0.3*	0–1	0

FRC: functional residual capacity; SD: standard deviation; BAL: bronchoalveolar lavage.

instillation of topical anesthesia in the segment for lavage, as this could inhibit bacterial growth.

Localization

In children with localized lung lesions, BAL should be performed in the most affected segment or lobe. In diffuse disease, any place may be chosen, although the most appropriate segments are the middle lobe (ML) and the lingula, both due to their accessibility and easy impaction of the bronchoscope as well as the greater recovery of the liquid instilled as it is favored by the pull of gravity (strong recommendation; high-quality evidence). In infants, one of the most accessible areas is the right lower lobe.

Liquid Used, Methods of Instillation and Recuperation

Sterile isotonic saline is used, at room temperature or at 37 °C to diminish the cough reflex.⁴⁴ It is introduced through the bronchoscope channel with a syringe, followed by the instillation of sufficient air in order to completely empty it. Afterwards, the liquid that was introduced is aspirated by either manual or mechanical aspiration through a circuit made up of several interconnected recipients between the working channel and the aspiration, applying an intermittent negative pressure between –25 and –100 mmH₂O (strong recommendation; low-quality evidence). Excessive negative pressure should be avoided as it could lead to the collapse of the distal airway, impeding the recovery of the liquid.

Liquid Volumes

Children <20 kg: a total of 3 ml/kg, divided into 3 equal fractions; children ≥20 kg: aliquots of 20 ml, up to a total final volume of 3 ml/kg (strong recommendation; low-quality evidence).

Processing of the Recovered Liquid

The volume instilled and the volume recovered both should be recorded. A recovery of 40% with few epithelial cells is accepted, although some authors⁴⁸ admit ≥10%. In healthy children, between 43% and 85% of the volume introduced can be recovered, although the first sample is usually clearly less than the remainder. The first aliquot is considered the bronchial sample. It has more neutrophils and less lymphocytes than the remainder. For this reason, it should be separated from the rest for microbiological cultures, unless it is necessary to analyze specifically the cell component in bronchial inflammatory processes like asthma. The other aliquots offer comparable results, and therefore can be joined in one container. The processing of the liquid obtained should be done within 4 h, and until processing it should be kept at 4 °C to favor cell viability.

Cellularity

Total cell count should be determined (considered valid with >300–500 cells/ml) as should the differential formula. Studies on the normal values for children do not include many cases, many of which are not always absolutely healthy, and they have different methodologies (Table 6).^{49–53}

Microorganisms

For the diagnosis of bacterial infections, Gram stain, quantitative colony-culture and, in some cases, intracellular bacteria studies should be done. The detection of more than 1% of epithelial squamous cells can indicate an excessive contamination of the sample by oropharyngeal flora. Isolations ≥10 000 UFC/ml are considered significant (strong recommendation; high-quality evidence). The detection of microorganisms that do not form part of the usual oropharyngeal flora (*M. tuberculosis*, *Legionella*, *Nocardia*) is always pathological (strong recommendation; high-quality evidence).

Diagnostic Indications

- Lung infections in immunocompromised children, those with transplantation or oncological diseases, who initiate with acute dyspnea and hypoxemia with new diffuse interstitial infiltrates on chest radiograph: in this case, the technique should be done before initiating antibiotic therapy or if there is no improvement after 48 h of treatment. The diagnostic performance is high at around 80%.
- Immunocompetent children with severe pneumonia who do not respond properly to treatment or who require mechanical ventilation, and in nosocomial pneumonia in intubated and ventilated children.
- Non-infectious lung diseases: this is the diagnostic method of choice for alveolar proteinosis, alveolar hemorrhage, pulmonary histiocytosis and pneumonia due to lipid accumulation. In the first, the milky appearance of BAL is characteristic, as are the presence of lipoprotein-like PAS (+) extracellular material, alcian blue (–) and the detection of the typical lamellar bodies by electron microscope. In the alveolar hemorrhage, the reddish appearance of the BAL, which increases in successive aliquots, induces suspicion. Perls stain identifies the hemosiderophages, although they may be absent in recent (less than 48 h) or remote (more than 12 days) hemorrhages. In adults, the stain is considered positive when there is more than 25% of macrophage stained, but in children a slight increase can already be significant. The identification of Langerhans cells induces suspicion of histiocytosis X. For its confirmation, the presence of $\geq 5\%$ of CD1a (+) cells or the positivity for S-100 protein is required.

Other diseases in which BAL findings are useful

- Acute eosinophilic pneumonia: eosinophils in the lavage $>25\%$.
- Pneumonitis due to hypersensitivity: there is lymphocytosis, as in adults, but with a normal CD4/CD8 relationship. It may be useful to evaluate the natural killer cells and the expression of antigen-DR in the human leukocytes.
- Aspiration syndromes: these can run their course with a high lipid index (percentage of macrophages loaded with lipids), although the threshold for positivity and its significance is up for debate.
- Interstitial lung disease: can provide information on the cell pattern and direct the diagnosis.
- Infants and toddlers with severe recurrent wheezing: the type of cells detected or certain inflammatory components (cysteinyl-leukotrienes) could lean towards a phenotype, asthmatic or not. At these ages, bronchial inflammation is usually basically neutrophilic, the same as in the severe bronchiolitis or bacterial bronchitis of children with chronic cough or recurring respiratory pathology.
- Cystic fibrosis: BAL in the first stages of the disease can reveal the existence of inflammation and/or bronchial infection.³¹

Therapeutic Indications

One accepted application is the treatment of alveolar proteinosis. In children under the age of 10–12, unlike in older children or in adults, it is not possible to use a double lumen endotracheal tube and ventilate one lung while another is washed. The standard technique involves placing, with the help of a flexible bronchoscope, a balloon catheter or a small endotracheal tube in one bronchus, inflate the balloon to isolate the bronchial tree and carry out a massive BAL through the tube or catheter, while ventilating the other lung by means of a laryngeal mask or an oropharyngeal tube (strong recommendation; low-quality evidence).

Bronchial Biopsy

This is done in an endobronchial lesion or in the mucosa. In this case, it allows for the epithelium, the basal membrane and, occasionally, the smooth muscle to be studied.⁵⁴ It can be used for the diagnosis of granulomatous diseases (tuberculosis or sarcoidosis), endobronchial tumors or chronic fungal diseases.⁵⁵ It is also appropriate for obtaining ciliated cells in the diagnosis of primary ciliary dyskinesia,⁵⁶ although nasal brushing, with equal performance and greater simplicity, is the most recommendable technique^{57,58} (strong recommendation; high-quality evidence); it is used in research studies in children with cystic fibrosis,^{59,60} asthma^{61–63} or recurrent wheezing bronchitis,^{64,65} and can be useful for the management of difficult-to-control asthma.⁶⁶

To carry out the procedure, standard 1.8 mm diameter forceps are used, utilizing the 4.9 mm bronchoscope with a 2–2.2 mm channel. The use of a laryngeal mask for the sedation and anesthesia allows for the use of this bronchoscope for biopsies in children over the ages of 2–2.5.⁶⁷ With the smaller 1.1 mm forceps, the sample is smaller, but some groups have achieved representative samples.⁶⁸ It is recommended to obtain 3–5 samples in the case of endobronchial lesions, and a minimum of 3, in the biopsies of the mucosa that are taken from a subsegmental bronchial carina, between the third and the fourth generation. The use of deep sedation or general anesthesia avoids patient movement during the procedure, increasing safety.⁶⁹

The technique is well tolerated and no important bleeding or pneumothorax has been reported. In biopsies of endobronchial lesions, the risk of bleeding increases, and if these are very vascularized, it is advisable to previously instill 1–2 ml of adrenalin at 1/20 000 and to carefully evaluate the risk–benefit relationship beforehand. Routine chest radiography is not necessary afterwards (strong recommendation; low-quality evidence).

Transbronchial Biopsy

With this technique, pulmonary parenchyma can be obtained for microscopic analysis, avoiding the need for thoracotomy. The majority of these procedures are carried out in lung transplant patients,⁷⁰ although it can also help diagnose other pediatric lung diseases.

This is the technique of choice for the diagnosis of rejection episodes in children with lung transplantation, with a sensitivity of 72%–94% and specificity of 90%–100% (strong recommendation; high-quality evidence). It has less value for the diagnosis of chronic rejection or bronchiolitis obliterans as in these cases the distribution of the lesions is very patchy and bronchioles are not always observed in the sample obtained. It should always be done when there is a suspicion for acute rejection, and selectively at 7–15 days, 6 weeks, 3, 6, and 12 months post-transplant.⁷¹ Some workgroups only indicate the procedure after 6 months in cases of suspected rejection.

It is not recommended for the diagnosis of idiopathic interstitial pneumopathies in children, in which case open biopsy is indicated⁷² (strong recommendation; high-quality evidence). It can be useful in granulomatous diseases, such as miliary tuberculosis or sarcoidosis, and also in extrinsic allergic alveolitis, eosinophilic pneumonia, vasculitis, alveolar proteinosis, histiocytosis X, alveolar microlithiasis, lymphocytic interstitial pneumonia and graft versus host disease.

In immunosuppressed patients,⁷³ transbronchial biopsy could help establish the definitive diagnosis in infections due to cytomegalovirus, by demonstrating the presence of intranuclear inclusions or a positive immunohistochemical stain in the lung tis-

sue, although given the focal nature of the infection, the frequency of positive results is only 11%–55%.

Before starting the test, it is necessary to have chest radiography, hemogram with platelets, coagulation tests, urea, creatinine and serum electrolytes. The procedure should be done with the patient under general anesthesia and either intubated or with a laryngeal mask, which gives better handling of the airway if the patient bleeds and makes for easier repetitive insertion and withdrawal of the bronchoscope. Afterwards, a follow-up radiograph is recommended within 2–4 h, in order to rule out the production of a pneumothorax and maintain the child in observation until the day after the procedure.

The use of fluoroscopy makes the technique easier, providing visualization and localization of the forceps, including their distance from the pleura, for the puncture⁷⁴ (strong recommendation; moderate-quality evidence).

With the 4.9 mm bronchoscope or with the rigid bronchoscope, standard 1.8 mm diameter adult biopsy forceps can be used. For the 3.6 mm or 2.8 mm bronchoscopes, 1.1 mm diameter smooth-edged biopsy forceps have been developed. These are used less often as the quantity of tissue that they extract is less and not as proper or standardized for histological studies. Nevertheless, it is possible to obtain representative samples in 85% of the procedures. Some authors prefer using the rigid bronchoscope for biopsies in small children as they accommodate the larger-sized standard forceps. However, it is more difficult to position the forceps in the segment desired and to access the upper lobes. The risk for pneumothorax is also greater.

If the affectation is focal, the biopsy is done in the altered lung lobe or segment. If the affectation is diffuse, it is usually done in the right lower lobe, which is the easiest and safest location. The middle lobe and lingula should be avoided due to the higher risk for pneumothorax. Biopsies will always be done in only one lung at each session, in order to avoid causing bilateral pneumothorax. It is reasonable to obtain 4–6 samples for the histological study of different segments of the biopsied lobe and one more sample if microbiological studies are necessary.

The most important complications are pneumothorax and hemorrhage. The incidence of pneumothorax oscillates between <1% and 3.5%²⁶ and it is more frequent in cases of bullous pulmonary lesions, peripheral localization of the lesions, mechanical ventilation and in immunosuppressed patients.⁷⁵ There is always a minimal amount of bleeding, which is self-limiting. Moderate (>25 ml) and severe (>100 ml) bleeding has been reported in 0.6%–5.4% of patients and, in some cases, death (1–2/1000 procedures).⁷³ It is important to properly select the patients and not perform biopsy on those with risk for bleeding: immunosuppression, uremia, lymphoma, leukemia, renal transplant and coagulation disorders.

Bronchial Brushing

This may be used to study the bronchial mucosa or to culture bronchial secretions. It is a useful technique for the diagnosis of ciliary dyskinesia, tuberculosis, pneumonia in patients with mechanical ventilation and lung infections in immunosuppressed patients.^{76,77}

Ideally, a protected brush should be used. It is made up of a brush housed in the interior of a catheter that, at the same time, is protected by another outer catheter. It is occluded in the distal end by a plug made up of an absorbable substance in order to prevent contamination on its way through the airway. The catheter is inserted through the channel of the bronchoscope until it surpasses the distal end by 3 cm in order to avoid collecting secretions accumulated in the tip of the catheter. The device is pushed to dis-

lodge the polyethylene glycol plug situated in at the distal end. The brush is moved outside the catheter and either the bronchial wall is brushed or the brush is rotated in the purulent secretions, if any are visualized, to take a sample. The brush is then retracted into the internal catheter, and then the internal catheter into the external catheter, and the bronchoscope is withdrawn.

The protected brush can only be used with the 4.9 mm bronchoscope. With the 3.6 or 2.8 mm bronchoscope, the non-protected brushes are used. In this case, upon withdrawing the brush through the working channel, the cells or the secretions may be lost. This risk diminishes if the brush is allowed to protrude from the bronchoscope (being careful not to touch the upper airway) or if the brush is hidden within the tip of the bronchoscope as it is being removed.

Once the brush is extracted, if it is protected, the outside is cleaned with 70% alcohol. If it is sent for culture, the brush is cut off with sterile scissors and it is deposited in a tube with 1 ml sterile physiological saline solution. If cytological or histological studies are ordered, the brush may be shaken in the appropriate medium to release the cells, spread on a slide or the tip cut and deposited in 10% formaldehyde.

Possible complications of bronchial brushing are bleeding (especially in patients with coagulation alterations or renal problems) and pneumothorax.

Rigid Bronchoscopy

Indications

FB has slowly been taking on a predominant role in pediatric bronchoscopy and, although it has substituted RB in most centers, rigid bronchoscopy is still widely used. The basic reason is the high incidence of foreign bodies in children. On the other hand, despite there being new ultrafine scopes for exploring newborns and infants, RB is still useful as a diagnostic or therapeutic tool when there is compromised ventilation, when extensive biopsies are necessary or when atelectasis should be resolved with the elimination of mucus plugging. The majority of the indications for RB are therapeutic⁷⁸: endoscopic treatment of localized airway obstruction, extraction of foreign bodies, management of massive hemoptysis and therapeutic instrumentation of the airway (strong recommendation; high-quality evidence). Pediatric RB is a procedure that is not regularly performed in most centers, therefore learning this technique usually requires specific training, generally postgraduate and in a reference center.

Obstruction of the Central Airways

The larynx, trachea and main bronchi are the bane of various surgical pathologies that cause stenosis of the lumen, from tumor to inflammatory lesions. RB is the endoscopic procedure of choice when open surgery cannot be contemplated to treat these processes either with the application of laser, implantation of endoprosthesis or other therapies. Although some of these techniques, such as the application of laser, balloon dilation or the placement of tracheobronchial prosthesis (except those made of silicone) can be done by bronchoscope or even by radioscopic control, there is extensive consensus that the safest method and that which gives the best results is rigid bronchoscopy.

In critical stenosis of the common airway with a reduction of $\leq 20\%$ of the predicted lumen and the patient's life at risk due to asphyxia, RB can be a life-saving procedure as it allows for the immediate restitution of the airway. It can dilate or perforate the nucleus of the tumor, opening the airway or progressively dilating the inflammatory stenosis. This not only avoids the need for tracheotomy if the obstruction is laryngeal, but it also the only option when the obstruction is located under the cervical trachea.

Extraction of Foreign Bodies

Rigid bronchoscopy is the procedure of choice for the safe and quick extraction of foreign bodies in children (strong recommendation; high-quality evidence). Nevertheless, in older children without compromised vital functions, FB, although it is more laborious and less effective (61% of successes compared with 97% for RB), can be an alternative.

Massive Hemoptysis

The majority of hemoptysis can be controlled with conservative measures, bronchoscopy and embolization by arteriography, but when the hemoptysis is massive and it compromises vital functions due to flooding of the airway with blood, RB is necessary for aspiration and for extracting the clots through its large working channel. Rigid bronchoscopy also accommodates two essential maneuvers at the same time: it guarantees ventilation while achieving hemostasis of the bleeding point. The former is done by protecting the contralateral lung from flooding by means of effective aspiration and selective intubation. Once the ventilation is guaranteed, the following objective is to achieve hemostasis. For this, either electrocoagulation or laser can be used. The hemorrhage can also be controlled by plugging the bleeding point with the bronchoscope itself or with an endobronchial occlusion balloon.

Therapeutic Instrumentation of the Airway

Among the therapies that are applied preferentially with RB are: laser, endoprosthesis, electrocoagulation, cryotherapy, argon plasma coagulation and balloon dilation. It can also be used as a guide in percutaneous tracheotomy and for the treatment of different lesions, such as tracheoesophageal fistulas or tracheo-bronchomalacia.

Instruments and Equipment

Types of Rigid Bronchoscopes

Pediatric bronchoscopes are shorter than those used in adults. They measure between 30 and 16 cm long with an internal diameter between 3.2 mm and 7 mm. To these measurements, one must also add the thickness of the wall of the tube, which ranges between 2 and 3 mm depending on the manufacturer, to obtain the external diameter. The beveled tip facilitates vision during intubation, and makes for easier passage through the vocal cords and stenosis and for mechanical resection. On its distal third, there are lateral holes, introduced by Killian, that ventilate the contralateral lung.

Some bronchoscope models have different ports of entry on their proximal end designed for different uses. Those of standard size are for connecting respirators and jet ventilation. Others allow for the simultaneous introduction of different types of fiber lasers and aspiration catheters. There are also optics of different thickness that adapt to the working channel.

In addition to the basic apparatus, it is convenient to have a video camera adaptable to the head of the optics. The variety of endoscopic forceps is extensive, but the main ones are: alligator mouth forceps for foreign bodies, universal forceps, spoon-shaped biopsy forceps and grasping forceps for soft foreign bodies.

Technique

The patient is positioned in supine decubitus with the bronchoscopist at his/her head. The ideal position for the visualization of the larynx and intubation is the combination of anterior flexion of the neck with hyperextension of the atlanto-occipital joint (flexion–extension or Boyde-Jackson position).

For intubation with the rigid bronchoscope, a laryngoscope can be used, using the standard tracheal intubation technique or directly using the RB; in this case, one may either look directly through the channel of the bronchoscope or use the straight optical viewfinder, which gives a better-quality magnified vision. The bronchoscope is held with the right hand and with the left the mouth is opened, protecting the incisors as well as the soft parts, lips and tongue from injury. The visualization of the arytenoids, with the vocal cords in front and the entrance of the esophagus behind, serves as a reference. In order to pass through the vocal cords, the bronchoscope is rotated on its axis clockwise until its tip coincides with the major axis of the *rima glottidis*. Once past the separated vocal cords, the rotation is undone. In the previously intubated patient, the point of the bronchoscope slips down on top of the orotracheal tube, using the tube as a guide until the tip comes to the point above the vocal cords; the orotracheal tube is then withdrawn and the RB pushed forward. Lastly, when there is a tracheotomy, the bronchoscope can also be introduced, carrying out what is known as a lower bronchoscopy in order to distinguish it from a translaryngeal or upper bronchoscopy. Once past the glottis, the ventilation system is connected to the lateral connections.

Complications

The complications of RB are due to the instrumentation with the bronchoscope itself, the medication used, the ventilation technique, the underlying pathology, the experience of the endoscopist and the type of intervention (extraction of foreign bodies, laser, placement of endoprosthesis, etc.).

The introduction of the bronchoscope can injure teeth and the soft parts of the mouth, the lips and tongue, as they can become caught between the metallic tube and the teeth. During the maneuver, the bronchoscopist must avoid using the incisors as a leverage point in order to elevate the soft tissue and expose the glottis. To prevent this type of injury, a dental protector is used. When the bronchoscope goes through the glottis, there may be damage to the arytenoids or vocal cords. When it enters into the trachea, there may be a cardiovascular response that, in adults, is usually sympathetic- tonic with an increase in the heart rate and blood pressure, but in younger people and in children the response may be vagotonic with bradycardia.

The most frequent complications are related with anesthesia and ventilation. These are hypoxemia, hypercapnia, cardiac arrhythmias and barotrauma.

Conflict of Interest

The authors declare having no conflict of interests.

References

- Schünemann HJ, Jaeschke R, Cook DJ, Bria WF, El-Solh AA, Ernst A, et al. An official ATS statement: grading the quality of evidence and strength of recommendations in ATS guidelines and recommendations. *Am J Respir Crit Care Med.* 2006;174:605–14.
- Wood RE. The emerging role of flexible bronchoscopy in pediatrics. *Clin Chest Med.* 2001;22:311–7.
- Pérez Ruiz E, Barrio Gómez de Agüero MI, Grupo de Técnicas de la Sociedad Española de Neumología Pediátrica. Broncoscopia flexible en el niño: indicaciones y aspectos generales. *An Pediatr.* 2004;60:354–66.
- Midulla F, De Blic J, Barbato A, Bush A, Eber E, Kotecha S, et al. Flexible endoscopy of paediatric airways. *Eur Respir J.* 2003;22:698–708.
- Barrio MI, Antelo C, Martínez MC. Fibrobroncoscopia. *An Pediatr Contin.* 2005;3:307–10.
- Nicolai T. Pediatric bronchoscopy. *Pediatr Pulmonol.* 2001;31:150–64.
- Wood RE. Evaluation of the upper airway in children. *Curr Opin Pediatr.* 2008;20:266–71.
- Figuerola Mulet J, Osona Rodríguez de Torres B, Llull Ferretjans M, Román Piñana JM. Contribución de la fibrobroncoscopia al diagnóstico de las enfermedades de la vía aérea superior. *An Pediatr.* 2005;63:137–42.

9. Boudewyns A, Claes J, Van de Heyning P. Clinical practice: an approach to stridor in infants and children. *Eur J Pediatr*. 2010;169:135–41.
10. Pérez Frías J, Caro Aguilera P, Pérez Ruiz E, Moreno Requena L. Tratamiento del cuerpo extraño intrabronquial. Broncoscopia combinada en neumología infantil. *An Pediatr (Barc)*. 2010;72:67–71.
11. Cakir E, Ersu RH, Uyan ZS, Oktem S, Karadag B, Yapar O, et al. Flexible bronchoscopy as a valuable tool in the evaluation of persistent wheezing in children. *Int J Pediatric Otorhinol*. 2009;73:1666–8.
12. Weinberger M, Abu-Hasan M. Pseudo-asthma: when cough, wheezing, and dyspnea are not asthma. *Pediatrics*. 2007;120:855–64.
13. Godfrey S. Pulmonary hemorrhage/hemoptysis in children. *Pediatr Pulmonol*. 2004;37:476–84.
14. Woods RK, Sharp RJ, Holcomb GW, Synder CL, Lofland GK, Ashcraft KW, et al. Vascular anomalies and tracheoesophageal compression: a single institutions 25-year experience. *Ann Thorac Surg*. 2001;72:434–8.
15. Efrati O, Gonik U, Modan-Moses D, Bielorai B, Barak A, Vardi A, et al. The role of flexible fiberoptic bronchoscopy in evaluation of pulmonary diseases in children with congenital cardiac disease. *Cardiol Young*. 2007;17:140–4.
16. Cerda J, Chacón J, Reichhard C, Bertrand P, Holmgren NL, Clavería D, et al. Flexible fiberoptic bronchoscopy in children with heart diseases: a twelve years experience. *Pediatr Pulmonol*. 2007;42:319–24.
17. Goussard P, Gie R. Airway involvement in pulmonary tuberculosis. *Pediatr Respir Rev*. 2007;8:118–23.
18. Cakir E, Uyan Z, Oktem S, Karakoc F, Ersu R, Karadag B, et al. Flexible bronchoscopy for diagnosis and follow up of childhood endobronchial tuberculosis. *Pediatr Infect Dis*. 2008;27:783–7.
19. Priftis KS, Anthracopoulos MB, Mermiri D, Papadopoulou A, Xepapadaki P, Tsakanika C, et al. Bronchial hyperresponsiveness, atopy and bronchoalveolar lavage eosinophils in persistent middle lobe syndrome. *Pediatr Pulmonol*. 2006;41:805–11.
20. Rubin BK. Respecting the medium lobe syndrome. *Pediatr Pulmonol*. 2006;41:803–4.
21. Priftis KS, Mermiri D, Papadopoulou A, Anthracopoulos MB, Vaos G, Polyxeni N. The role of timely intervention in middle lobe syndrome in children. *Chest*. 2005;128:2504–10.
22. Masters IB, Zimmerman PV, Chang AB. Longitudinal Quantification of growth and changes in primary tracheobronchomalacia sites in children. *Pediatr Pulmonol*. 2007;42:906–13.
23. Wright CD. Mediastinal tumors and cysts in the pediatric population. *Thorac Surg Clin*. 2009;19:47–61.
24. Nussbaum E. Pediatric fiberoptic bronchoscopy: clinical experience with 2.836 bronchoscopies. *Pediatr Crit Care Med*. 2002;3:171–6.
25. Eber E, Oberwaldner B. Tracheostomy care in the hospital. *Paediatr Respir Rev*. 2006;7:175–84.
26. McShane D, Nicholson AG, Goldstraw P, Ladas G, Travis WD, Ramanan R, et al. Inflammatory endobronchial polyps in childhood: clinical spectrum and possible link to mechanical ventilation. *Pediatr Pulmonol*. 2002;34:79–84.
27. Pérez-Ruiz E, Pérez-Frías J, Martínez-González B, Martínez-Arán T. Fibrobroncoscopia pediátrica. Análisis de una década. *An Esp Pediatr*. 2001;55:421–8.
28. Pérez Frías J, Pérez Ruiz E, Cordón Martínez A. Pulmonary diagnostic procedures: bronchoscopy. *Pediatr Pulmonol*. 2004;26S:234–7.
29. Saito J, Harris WT, Gelfonf J, Noah TL, Leigh MW, Johnson R, et al. Physiologic, bronchoscopic and bronchoalveolar lavage fluid findings in young children with recurrent wheeze and cough. *Pediatr Pulmonol*. 2006;41:709–19.
30. Barret JP, Desai MH, Hemdon DN. Effects of tracheostomies on infection and airway complications in pediatric burn patients. *Burns*. 2000;26:190–3.
31. Barrio Gómez de Agüero MI. Fibrosis quística: fibrobroncoscopia y BAL en la enfermedad broncopulmonar. *An Pediatr (Barc)*. 2009;71S:24–6.
32. Dabo L, Qiyi Z, Jianwen Z, Zhenyun H, Lifeng Z. Perioperative management of plastic bronchitis in children. *Inter J Pediatr Otorrhinol*. 2010;74:15–21.
33. Riethmueller J, Borth-Bruhns T, Kumpf M, Vonthein R, Wiskirchen J, Stern M, et al. Recombinant human deoxyribonuclease shortens ventilation time in young, mechanically ventilated children. *Pediatr Pulmonol*. 2006;41:61–6.
34. Xu Z, Jing J, Wang H, Xu F, Wang J. Pulmonary alveolar proteinosis in China: a systematic review of 241 cases. *Respirology*. 2009;14:761–6.
35. Fiadjo J, Stricker P. Pediatric difficult airway management: current devices and techniques. *Anesthesiol Clin*. 2009;27:185–95.
36. Tang LF, Xu YC, Wang YS, Wang CF, Zhu GH, Bao XE, et al. Airway foreign body removal by flexible bronchoscopy: experience with 1027 children during 2000–2008. *World J Pediatr*. 2009;5:191–5.
37. Lain A, Fanjul M, García Casillas MA, Parente A, Cañizo A, Carreras N, et al. Extracción de cuerpos extraños en vía aérea pediátrica con fibrobroncoscopia. *Cir Pediatr*. 2007;20:194–8.
38. Rignini CA, Morel N, Karkas A, Reyt E, Ferreti K, Pin I, et al. What is the diagnosis value of flexible bronchoscopy in the initial investigation of children with suspected foreign body aspiration? *Int J Pediatr Otorhinolaryngol*. 2007;71:1383–90.
39. Meier JD, Sulman CG, Almond PS, Holinger LD. Endoscopic management of recurrent congenital tracheoesophageal fistula: a review of techniques and results. *Int J Pediatr Otorrhinol*. 2007;71:691–7.
40. Atzori P, Iacobelli BD, Bottero S, Spiriyidakis J, Laviani R, Trucchi A, et al. Pre-operative tracheobronchoscopy in newborns with esophageal atresia: does it matter? *J Pediatr Surg*. 2006;41:1054–7.
41. Fayon M, Donato L, De Blic J, Labbé A, Becmeur F, Mely L, et al. French experience of silicone tracheobronchial stenting in children. *Pediatr Pulmonol*. 2005;39:21–7.
42. De Blic J, Marchac V, Scheinmann P. Complications of flexible bronchoscopy in children: prospective study of 1.328 procedures. *Eur Respir J*. 2002;20:1271–6.
43. Solomonov A, Fruchter O, Zuckerman T, Brenner B, Yigla M. Pulmonary hemorrhage: a novel mode of therapy. *Respir Med*. 2009;103:1196–200.
44. Escribano Montaner A, Moreno Galdó A, Grupo de Técnicas de la Sociedad Española de Neumología Pediátrica. Técnicas fibrobroncoscópicas especiales: lavado broncoalveolar, biopsia bronquial y biopsia transbronquial. *An Pediatr (Barc)*. 2005;62:352–66.
45. De Blic J, Delacourt C, Scheinmann P. Ultrathin flexible bronchoscopy in neonatal intensive care unit. *Arch Dis Child*. 1991;66:1383–5.
46. Hasegawa S, Hitomi S, Murakawa M, Mori K. Development of an ultrathin fiberscope with a built-in channel for bronchoscopy in infants. *Chest*. 1996;110:1543–6.
47. Pérez Frías J, Pérez Ruiz E, Cordón Martínez A, Spitaleri G, editors. *Broncoscopia pediátrica*. 2nd ed. Madrid: Ergon; 2008. p. 65–83.
48. Efrati O, Sadeh-Gornik U, Modan-Moses D, Barak A, Szeinberg A, Vardi A, et al. Flexible bronchoscopy and bronchoalveolar lavage in pediatric patients with lung disease. *Pediatr Crit Care Med*. 2009;10:80–4.
49. Clement A, Chadelat K, Masliah J, Housset B, Sartet A, Grimfeld A, et al. A controlled study of oxygen metabolite release by alveolar macrophages from children with interstitial lung disease. *Am Rev Respir Dis*. 1987;136:1424–8.
50. Ratjen F, Bredendiek M, Brendel M, Meltzer J, Costabel U. Differential cytology of bronchoalveolar lavage fluid in normal children. *Eur Respir J*. 1994;7:1865–70.
51. Riedler J, Grigg J, Stone C, Tauro G, Robertson CF. Bronchoalveolar lavage cellularity in healthy children. *Am J Respir Crit Care Med*. 1995;152:163–8.
52. Midulla F, Villani A, Merolla R, Bjermer L, Sandstrom T, Ronchetti R. Bronchoalveolar lavage studies in children without parenchymal lung disease: cellular constituents and protein levels. *Pediatr Pulmonol*. 1995;20:112–8.
53. Tessier V, Chadelat K, Baculard A, Housset B, Clement A. A controlled study of differential cytology and cytokine expression profiles by alveolar cells in pediatric sarcoidosis. *Chest*. 1996;109:1430–8.
54. Bush A, Pohunek P. Brush biopsy and mucosal biopsy. *Am J Respir Crit Care Med*. 2000;162:S18–22.
55. Whitehead BF. Bronchoscopic lung biopsy in pediatric patients. *J Bronchol*. 1999;6:48–54.
56. Barbato A, Frischer T, Kuehni CE, Snijders D, Azevedo I, Baktai G, et al. Primary ciliary dyskinesia: a consensus statement on diagnostic and treatment approaches in children. *Eur Respir J*. 2009;34:1264–76.
57. MacCormick J, Robb I, Kovess T, Carpenter B. Optimal biopsy techniques in the diagnosis of primary ciliary dyskinesia. *J Otolaryngol*. 2002;31:13–7.
58. Nuesslein TG, Hufnagel C, Stephan V, Rieger CH. Yield of bronchial biopsies in addition to nasal brushing for ciliary function analyses in children. *Klin Padiatr*. 2004;216:238–43.
59. Regamey N, Hilliard TN, Saglani S, Zhu J. Quality, size, and composition of pediatric endobronchial biopsies in cystic fibrosis. *Chest*. 2007;131:1710–7.
60. Molina-Teran A, Hilliard TN, Saglani S, Haxby E. Safety of endobronchial biopsy in children with cystic fibrosis. *Pediatr Pulmonol*. 2006;41:1021–4.
61. Cokugras H, Akcakaya N, Seckin I, Camcioglu Y, Sarimurat N, Aksoy F. Ultrastructural examination of bronchial biopsy specimens from children with moderate asthma. *Thorax*. 2001;56:25–9.
62. Payne DNR, Adcock IM, Wilson NM, Oates T, Scallan M, Bush A. Relationship between exhaled nitric oxide and mucosal eosinophilic inflammation in children with difficult asthma, after treatment with oral prednisolone. *Am J Respir Crit Care Med*. 2001;164:1376–81.
63. Payne DNR, Rogers AV, Adelroth E, Vandi V, Guntupalli KK, Bush A, et al. Early thickening of the reticular basement membrane in children with difficult asthma. *Am J Respir Crit Care Med*. 2003;167:78–82.
64. Elston WJ, Whittaker AJ, Khan LN, Flood-Page P. Safety of research bronchoscopy, biopsy and bronchoalveolar lavage in asthma. *Eur Respir J*. 2004;24:375–7.
65. Saglani S, Nicholson AG, Scallan M, Balfour-Lynn I, Rosenthal M, Payne DN, et al. Investigation of young children with severe recurrent wheeze: any clinical benefit? *Eur Respir J*. 2006;27:29–35.
66. De Blic J, Tillie-Leblond I, Tonnel AB, Jaubert F, Scheinmann P, Gosset P. Difficult asthma in children: an analysis of airway inflammation. *J Allergy Clin Immunol*. 2004;113:94–100.
67. Salva PS, Theroux C, Schwartz D. Safety of endobronchial biopsy in 170 children with chronic respiratory symptoms. *Thorax*. 2003;58:1058–60.
68. Saglani S, Payne DN, Nicholson AG, Scallan M. The safety and quality of endobronchial biopsy in children under five years old. *Thorax*. 2003;58:1053–7.
69. Payne DN, McKenzie SA, Stacey S, Misra D, Haxby E, Bush A. Safety and ethics of bronchoscopy and endobronchial biopsy in difficult asthma. *Arch Dis Child*. 2001;84:423–6.
70. Visner GA, Faro A, Zander DS. Role of transbronchial biopsies in pediatric lung diseases. *Chest*. 2004;126:273–80.
71. Benden B, Harpur-Sinclair O, Ranasinghe AS, Hartley JC. Surveillance bronchoscopy in children during the first year after lung transplantation: is it worth it? *Thorax*. 2007;62:57–61.
72. Fan LL, Kozinetz CA, Wojczak HA, Chatfield BA, Cohen AH, Rothenberg SS. Diagnostic value of transbronchial, thorascopic, and open lung biopsy in immunocompetent children with chronic interstitial lung disease. *J Pediatr*. 1997;131:565–9.
73. Becker HD, Shirakawa T, Tanaka F, Müller KM, Herth F. Transbronchial lung biopsy in the immunocompromised patient. *Eur Respir Mon*. 1998;9:193–208.

74. Whitehead B, Scott JP, Helms P, Malone M, Macrae D, Higenbottam TW, et al. Technique and use of transbronchial biopsy in children and adolescents. *Pediatr Pulmonol.* 1992;12:240–6.
75. McDougall JC, Cortese DA. Bronchoscopic lung biopsy. In: Prakash UBS, editor. *Bronchoscopy.* New York: Raven Press Ltd; 1994. p. 141–6.
76. Ioanas M, Ferrer R, Angrill J, Ferrer M, Torres A. Microbial investigation in ventilator-associated pneumonia. *Eur Respir J.* 2001;17:791–801.
77. Rañó A, Agustí C, Jiménez P, Angrill J, Benito N, Danés C, et al. Pulmonary infiltrates in non-HIV immunocompromised patients: a diagnostic approach using non-invasive and bronchoscopic procedures. *Thorax.* 2001;56:379–87.
78. Díaz-Agero P, Canseco F, Gil JL. Indicaciones y técnica de la broncoscopia rígida. In: Díaz-Agero P, Flandes J, editors. *Broncoscopia diagnóstica y terapéutica.* Madrid: Ergon; 2007. p. 53–70.